

[Open Peer Review on Qeios](#)

# Luxol Fast Blue and Cresyl Violet Staining Method

National Cancer Institute

## Source

National Cancer Institute. *Luxol Fast Blue and Cresyl Violet Staining Method*. NCI Thesaurus. Code C92300.

A microscopy staining method that utilizes luxol fast blue, an alcohol soluble copper phthalocyanine dye, and cresyl violet, an organic compound that is commonly used as a counter stain, to determine the presence of myelin as well as helicobacter and campylobacter bacterial species.